

# Computational Methods for Biomedical Image Analysis

Assignment 1

Medical Image Segmentation


# Introduction

- In this assignment, we will learn
  - Usage of DICOM (Digital Imaging and Communications in Medicine) files
  - Hounsfield Units
  - CT Segmentation using thresholding algorithm

# Part 1 (10%)

- Read in and print out all the data fields in a DICOM file (one slice) (5%)
- Read in the raw data for a CT slice and convert its pixel values into Hounsfield units. Compute the max, min, mean and standard deviation of both images (raw data and Hounsfield units) (5%)

Hounsfield Units =  $RescaleSlope \times Raw\ Image + RescaleIntercept$




Bone	+400 → +1000
Soft tissue	+40 → +80
Water	0
Fat	-60 → -100
Lung	-400 → -600
Air	-1000

```
(0008, 0005) Specific Character Set      CS: 'ISO_IR 100'
(0008, 0016) SOP Class UID               UI: CT Image Storage
(0008, 0018) SOP Instance UID            UI: 1.3.6.1.4.1.14519.5.2.1.7009.9004.993926557618716725548364243182
(0008, 0060) Modality                    CS: 'CT'
(0008, 103e) Series Description           LO: 'Axial'
(0010, 0010) Patient's Name               PN: '0acbebb8d463b4b9ca88cf38431aac69'
(0010, 0020) Patient ID                  LO: '0acbebb8d463b4b9ca88cf38431aac69'
(0010, 0030) Patient's Birth Date        DA: '19000101'
(0018, 0060) KVP                         DS: ''
(0020, 000d) Study Instance UID          UI: 2.25.66234994940093060530875882673593880723182397297085825139120
(0020, 000e) Series Instance UID         UI: 2.25.27985737130106072918310533525688877208529713445697698517643
(0020, 0011) Series Number               IS: '4'
(0020, 0012) Acquisition Number          IS: '2'
(0020, 0013) Instance Number             IS: '1'
(0020, 0032) Image Position (Patient)    DS: ['-157.2275390625', '-280.2275390625', '-18.9']
(0020, 0037) Image Orientation (Patient) DS: ['1', '0', '0', '0', '1', '0']
(0020, 0052) Frame of Reference UID      UI: 2.25.11024047567322936673371602952341875417234268835153580015625
(0020, 1040) Position Reference Indicator LO: ''
(0020, 1041) Slice Location              DS: "-18.9"
(0028, 0002) Samples per Pixel            US: 1
(0028, 0004) Photometric Interpretation  CS: 'MONOCHROME2'
(0028, 0010) Rows                        US: 512
(0028, 0011) Columns                     US: 512
(0028, 0030) Pixel Spacing                DS: ['0.544921875', '0.544921875']
(0028, 0100) Bits Allocated               US: 16
(0028, 0101) Bits Stored                  US: 12
(0028, 0102) High Bit                     US: 11
(0028, 0103) Pixel Representation         US: 0
(0028, 0106) Smallest Image Pixel Value   US: 0
(0028, 0107) Largest Image Pixel Value    US: 4090
(0028, 0301) Burned In Annotation         CS: 'NO'
(0028, 0303) Longitudinal Temporal Information M SH: 'MODIFIED'
(0028, 1050) Window Center                DS: ['-400', '0']
(0028, 1051) Window Width                 DS: ['1600', '2']
(0028, 1052) Rescale Intercept            DS: "-1024"
(0028, 1053) Rescale Slope                DS: "1"
(0028, 1055) Window Center & Width Explanation LO: ['WINDOW1', 'WINDOW2']
(7fe0, 0010) Pixel Data                   OW: Array of 524288 bytes
```

## Part 2 (20%)

- In part 1, we read in only one individual slice (a DICOM file). Now we want to read in a 3D volume. (DICOM files in a folder)
- Note that you will need to sort all the slices to make it into correct order. Please explain how do you sort the slices. (10%)
- Normalize all the pixel from Hounsfield Units to float32 type number between 0.0 to 1.0 and display 25 slices in correct order. (10%)

## Part 2 (Cont.)



## Part 3 (60%)




- Now we can try to segment out our lung. We know that lower Hounsfield units correspond to low density materials (like air) and higher Hounsfield units correspond to highly attenuative materials, like bone.
- Please try to use at least two different thresholding algorithm to segment the lung.
- Threshold Value
  - Balanced Histogram Thresholding (BHT)
  - Local mean
  - Local median
  - Otsu's method

### Thresholding Algorithm

- Choose a threshold pixel value  $T$
- For every pixel
  - if pixel  $\geq T$ , label as foreground
  - else label as background


# Part 3 (Cont.)

- Each threshold method: (30%)
  - Plot the histogram of your pixels and the threshold. (15%)
  - Display one CT slice and the corresponded segmentation result. (15%)



# Bonus (10%)

- Try to show your best 3D segmentation result with 3D Viewer.





# Environment Requirement and Dataset

- Python 3
- Python package:
  - Numpy
  - Pydicom
  - Matplotlib
  - Skimage
- Please download the dataset from the following link:  
[https://drive.google.com/file/d/1AcO7rZIpSJhPLjgl5ghLxJ0uzqdbBI\\_w/view?usp=sharing](https://drive.google.com/file/d/1AcO7rZIpSJhPLjgl5ghLxJ0uzqdbBI_w/view?usp=sharing)

# Report (10%)

- Please write your report with the given latex template in English.
- Show your result inside the report and give some explanations.
- Please also give a short summary about what you have learned.

# Reminder

- Please zip your code (.py) and report (.pdf) together and name as {Student\_id}\_{Name}.zip
- **Deadline: Upload to iLMS before 2020/3/29 23:59 (UTC+8).**
- Late submission: 20% penalty.
- Wrong hand in format: 10% penalty.
- **This is an individual assignment, no cheating.**

# Hint

- You can write your report through Overleaf, an online latex editor.
- You can write your code with jupyter notebook.